

WDFY1 Polyclonal Antibody

Cat No: HR1AP12854

For research use only

Overview

Product Name	WDFY1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	45kDa
GeneID?Human?	WDFY1 KIAA1435 WDF1 ZFYVE17
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	WD repeat and FYVE domain containing 1(WDFY1) Homo sapiens The protein encoded by this gene is a phosphatidylinositol 3-phosphate binding protein, which contains a FYVE zinc finger domain and multiple WD-40 repeat domains. When exogenously expressed, it localizes to early endosomes. Mutagenesis analysis demonstrates that this endosomal localization is mediated by the FYVE domain. [provided by RefSeq, Jan 2015],