

GPC39 Polyclonal Antibody

Cat No: HR1AP12845

For research use only

Overview

Product Name	GPC39 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	29kDa
GeneID?Human?	
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	GPR32P1 (G Protein-Coupled Receptor 32, Pseudogene 1) is a Pseudogene. GO annotations related to this gene include G-protein coupled receptor activity.