

TAAR8 Polyclonal Antibody

Cat No: HR1AP12843

For research use only

Overview

Product Name	TAAR8 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	TAAR8 GPR102 TA5 TAR5 TRAR5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	trace amine associated receptor 8(TAAR8) Homo sapiens This gene is part of the trace amine receptor cluster on chromosome 6 and encodes an orphan G-protein coupled receptor. Upregulated expression of this gene in astroglial cells upon exposure to lipopolysaccharides suggests a function for the encoded protein in the brain. [provided by RefSeq, Jul 2016],