

NPFF2 Polyclonal Antibody

Cat No: HR1AP12770

For research use only

Overview

Product Name	NPFF2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	57kDa
GeneID?Human?	NPFFR2 GPR74 NPFF2 NPGPR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	neuropeptide FF receptor 2(NPFFR2) Homo sapiens This gene encodes a member of a subfamily of G-protein-coupled neuropeptide receptors. This protein is activated by the neuropeptides A-18-amide (NPAF) and F-8-amide (NPFF) and may function in pain modulation and regulation of the opioid system. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2009],