



## MPRB Polyclonal Antibody

Cat No: HR1AP12767

For research use only

### Overview

Product Name	MPRB Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	PAQR8 C6orf33 LMPB1 MPRB
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Steroid membrane receptor. Binds progesterone. May be involved in oocyte maturation.,similarity:Belongs to the ADIPOR family.,subcellular location:Colocalizes with a lysosomal protein cathepsin D.,tissue specificity:Expressed in brain and spinal cord.,