

LT4R2 Polyclonal Antibody

Cat No: HR1AP12765

For research use only

Overview

Product Name	LT4R2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	LTB4R2 BLT2R BLTR2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>caution:It is uncertain whether Met-1 or Met-32 is the initiator.,function:Low-affinity receptor for leukotrienes including leukotriene B4. Mediates chemotaxis of granulocytes and macrophages. The response is mediated via G-proteins that activate a phosphatidylinositol-calcium second messenger system. The rank order of affinities for the leukotrienes is LTB4 > 12-epi-LTB4 > LTB5 > LTB3.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Widely expressed.,</p>