

VN1R1 Polyclonal Antibody

Cat No: HR1AP12762

For research use only

Overview

Product Name	VN1R1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	VN1R1 V1RL1 VNR1911
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	vomer nasal 1 receptor 1(VN1R1) Homo sapiens Pheromones are chemical signals that elicit specific behavioral responses and physiologic alterations in recipients of the same species. The protein encoded by this gene is similar to pheromone receptors and is primarily localized to the olfactory mucosa. An alternate splice variant of this gene is thought to exist, but its full length nature has not been determined. [provided by RefSeq, Jul 2008],