

## UR2R Polyclonal Antibody

Cat No: HR1AP12761

For research use only

### Overview

Product Name	UR2R Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	UTS2R GPR14
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:High affinity receptor for urotensin-2 and urotensin-2B. The activity of this receptor is mediated by a G-protein that activate a phosphatidylinositol-calcium second messenger system.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Most abundant expression in the heart and pancreas.,