

GP111 Polyclonal Antibody

Cat No: HR1AP12764

For research use only

Overview

| | |
|-----------------------|---|
| Product Name | GP111 Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 77kDa |
| GeneID?Human? | GPR111 PGR20 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | function:Orphan receptor.,similarity:Belongs to the G-protein coupled receptor 2 family. LN-TM7 subfamily.,similarity:Contains 1 GPS domain., |