

TSHR Polyclonal Antibody

Cat No: HR1AP12760

For research use only

Overview

Product Name	TSHR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	84kDa
GeneID?Human?	TSHR LGR3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>thyroid stimulating hormone receptor(TSHR) Homo sapiens The protein encoded by this gene is a membrane protein and a major controller of thyroid cell metabolism. The encoded protein is a receptor for thyrothopin and thyrostimulin, and its activity is mediated by adenylate cyclase. Defects in this gene are a cause of several types of hyperthyroidism. Three transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Dec 2008],</p>