

TS1R2 Polyclonal Antibody

Cat No: HR1AP12759

For research use only

Overview

Product Name	TS1R2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	92kDa
GeneID?Human?	TAS1R2 GPR71 T1R2 TR2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Putative taste receptor. TAS1R2/TAS1R3 recognizes diverse natural and synthetic sweeteners.,online information:The taste experience -Issue 55 of February 2005,similarity:Belongs to the G-protein coupled receptor 3 family. TAS1R subfamily.,subunit:Forms heterodimers with TAS1R3.,