

## TM7S3 Polyclonal Antibody

Cat No: HR1AP12757

For research use only

### Overview

Product Name	TM7S3 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	62kDa
GeneID?Human?	TM7SF3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	tissue specificity:Ubiquitous. The highest expression is in kidney.,