

TAAR6 Polyclonal Antibody

Cat No: HR1AP12754

For research use only

Overview

Product Name	TAAR6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	TAAR6 TA4 TAR4 TRAR4
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	trace amine associated receptor 6(TAAR6) Homo sapiens This gene encodes a seven-transmembrane G-protein-coupled receptor that likely functions as a receptor for endogenous trace amines. Mutations in this gene may be associated with schizophrenia.[provided by RefSeq, Feb 2010],