

# TAAR9 Polyclonal Antibody

Cat No: HR1AP12755

For research use only

## Overview

Product Name	TAAR9 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	TAAR9 TA3 TAR3 TRAR3
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	trace amine associated receptor 9 (gene/pseudogene)(TAAR9) Homo sapiens TAAR9 is a member of a large family of rhodopsin G protein-coupled receptors (GPCRs, or GPRs). GPCRs contain 7 transmembrane domains and transduce extracellular signals through heterotrimeric G proteins.[supplied by OMIM, Jul 2005],