

RGDSR Polyclonal Antibody

Cat No: HR1AP12743

For research use only

Overview

Product Name	RGDSR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	52kDa
GeneID?Human?	RGL4 RGR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>ral guanine nucleotide dissociation stimulator like 4(RGL4) Homo sapiens This oncogene encodes a protein similar to guanine nucleotide exchange factor Ral guanine dissociation stimulator. Increased expression of this gene leads to translocation of the encoded protein to the cell membrane. The encoded protein can activate several pathways, including the Ras-Raf-MEK-ERK cascade. [provided by RefSeq, Jul 2016],</p>