

PKR2 Polyclonal Antibody

Cat No: HR1AP12740

For research use only

Overview

Product Name	PKR2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	42kDa
GeneID?Human?	PROKR2 GPR73L1 PKR2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>prokineticin receptor 2(PROKR2) Homo sapiens Prokineticins are secreted proteins that can promote angiogenesis and induce strong gastrointestinal smooth muscle contraction. The protein encoded by this gene is an integral membrane protein and G protein-coupled receptor for prokineticins. The encoded protein is similar in sequence to GPR73, another G protein-coupled receptor for prokineticins. [provided by RefSeq, Jul 2008],</p>