

# PSYR Polyclonal Antibody

Cat No: HR1AP12741

For research use only

## Overview

Product Name	PSYR Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	GPR65 TDAG8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>function:Receptor for the glycosphingolipid psychosine (PSY) and several related glycosphingolipids. May have a role in activation-induced cell death or differentiation of T-cells.,similarity:Belongs to the G-protein coupled receptor 1 family.,tissue specificity:Predominantly expressed in thymus, spleen, lymph nodes, small intestine, lung, placenta and peripheral blood leukocytes.,</p>