

PACR Polyclonal Antibody

Cat No: HR1AP12737

For research use only

Overview

| | |
|-----------------------|--|
| Product Name | PACR Polyclonal Antibody |
| Source | Rabbit |
| Applications | WB,ELISA |
| Species Reactivity | Human,Rat,Mouse |
| Recommended Dilutions | |
| Immunogen | |
| Species | Rabbit |
| Storage | -20°C/1 year |
| Isotype | |
| Clonality | |
| Concentration | 1 mg/ml |
| Observed band | 51kDa |
| GeneID?Human? | ADCYAP1R1 |
| Human Swiss-Prot No. | |
| Cellular localization | |
| Alternative Names | |
| Background | ADCYAP receptor type I(ADCYAP1R1) Homo sapiens This gene encodes type I adenylate cyclase activating polypeptide receptor, which is a membrane-associated protein and shares significant homology with members of the glucagon/secretin receptor family. This receptor mediates diverse biological actions of adenylate cyclase activating polypeptide 1 and is positively coupled to adenylate cyclase. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified. [provided by RefSeq, Dec 2010], |