

P2Y10 Polyclonal Antibody

Cat No: HR1AP12735

For research use only

Overview

Product Name	P2Y10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	P2RY10
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>purinergic receptor P2Y10(P2RY10) Homo sapiens The protein encoded by this gene belongs to the family of G-protein coupled receptors that are preferentially activated by adenosine and uridine nucleotides. There is a pseudogene for this gene nearby on chromosome X. Multiple alternatively spliced transcripts have been observed. [provided by RefSeq, Apr 2016],</p>