

# LPAR4 Polyclonal Antibody

Cat No: HR1AP12734

For research use only

## Overview

Product Name	LPAR4 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	40kDa
GeneID?Human?	LPAR4 GPR23 LPA4 P2RY9
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lysophosphatidic acid receptor 4(LPAP4) Homo sapiens This gene encodes a member of the lysophosphatidic acid receptor family. It may also be related to the P2Y receptors, a family of receptors that bind purine and pyrimidine nucleotides and are coupled to G proteins. The encoded protein may play a role in monocytic differentiation. [provided by RefSeq, Feb 2009],