

LPAR6 Polyclonal Antibody

Cat No: HR1AP12732

For research use only

Overview

Product Name	LPAR6 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	LPAR6 P2RY5
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	lysophosphatidic acid receptor 6(LPAR6) Homo sapiens The protein encoded by this gene belongs to the family of G-protein coupled receptors, that are preferentially activated by adenosine and uridine nucleotides. This gene aligns with an internal intron of the retinoblastoma susceptibility gene in the reverse orientation. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jun 2009],