

MCHR1 Polyclonal Antibody

Cat No: HR1AP12618

For research use only

Overview

Product Name	MCHR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	46kDa
GeneID?Human?	MCHR1 GPR24 SLC1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>melanin concentrating hormone receptor 1(MCHR1) Homo sapiens The protein encoded by this gene, a member of the G protein-coupled receptor family 1, is an integral plasma membrane protein which binds melanin-concentrating hormone. The encoded protein can inhibit cAMP accumulation and stimulate intracellular calcium flux, and is probably involved in the neuronal regulation of food consumption. Although structurally similar to somatostatin receptors, this protein does not seem to bind somatostatin. [provided by RefSeq, Jul 2008].</p>