

GPR98 Polyclonal Antibody

Cat No: HR1AP12613

For research use only

Overview

Product Name	GPR98 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	693kDa
GeneID?Human?	GPR98 KIAA0686 KIAA1943 MASS1 VLGR1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>adhesion G protein-coupled receptor V1(ADGRV1) Homo sapiens This gene encodes a member of the G-protein coupled receptor superfamily. The encoded protein contains a 7-transmembrane receptor domain, binds calcium and is expressed in the central nervous system. Mutations in this gene are associated with Usher syndrome 2 and familial febrile seizures. Several alternatively spliced transcripts have been described. [provided by RefSeq, Jul 2008],</p>