

HCAR1 Polyclonal Antibody

Cat No: HR1AP12611

For research use only

Overview

Product Name	HCAR1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	38kDa
GeneID?Human?	HCAR1 GPR104 GPR81 HCA1 FKSG80
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	hydroxycarboxylic acid receptor 1(HCAR1) Homo sapiens G protein-coupled receptors (GPCRs, or GPRs), such as GPR81, contain 7 transmembrane domains and transduce extracellular signals through heterotrimeric G proteins.[supplied by OMIM, Feb 2005],