

GPR64 Polyclonal Antibody

Cat No: HR1AP12610

For research use only

Overview

Product Name	GPR64 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	111kDa
GeneID?Human?	GPR64 HE6 TM7LN2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>adhesion G protein-coupled receptor G2(ADGRG2) Homo sapiens This gene encodes a member of the G protein-coupled receptor family described as an epididymis-specific transmembrane protein. The encoded protein may be proteolytically processed as it contains a motif shown to be a protein scission motif in some members of this family (PMID: 11973329). Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Aug 2011],</p>