

GP125 Polyclonal Antibody

Cat No: HR1AP12605

For research use only

Overview

Product Name	GP125 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	145kDa
GeneID?Human?	GPR125 UNQ556/PRO1113
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>adhesion G protein-coupled receptor A3(ADGRA3) Homo sapiens This gene encodes a member of the G protein-coupled receptor superfamily. This membrane protein may play a role in tumor angiogenesis through its interaction with the human homolog of the Drosophila disc large tumor suppressor gene. This gene is mapped to a candidate region of chromosome 4 which may be associated with bipolar disorder and schizophrenia. [provided by RefSeq, Oct 2012].</p>