

GP141 Polyclonal Antibody

Cat No: HR1AP12606

For research use only

Overview

Product Name	GP141 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	GPR141 PGR13
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G protein-coupled receptor 141(GPR141) Homo sapiens GPR141 is a member of the rhodopsin family of G protein-coupled receptors (GPRs) (Fredriksson et al., 2003 [PubMed 14623098]).[supplied by OMIM, Mar 2008],