

## GP112 Polyclonal Antibody

Cat No: HR1AP12604

For research use only

### Overview

Product Name	GP112 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	338kDa
GeneID?Human?	GPR112
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>adhesion G protein-coupled receptor G4(ADGRG4) Homo sapiens This gene encodes a G-protein coupled receptor belonging to a large family of diverse integral membrane proteins that participate in various physiological functions. Members of this superfamily are characterized by a signature 7-transmembrane domain motif. The ligand for this family member is unknown, and it is therefore an orphan receptor. This receptor is known to be expressed in normal enterochromaffin cells and in gastrointestinal neuroendocrine carcinoma cells, and it is therefore considered to be a novel biomarker or target for immunotherapy. [provided by RefSeq, May 2010],</p>