

## GLP2R Polyclonal Antibody

Cat No: HR1AP12602

For research use only

### Overview

Product Name	GLP2R Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	GLP2R
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	glucagon like peptide 2 receptor(GLP2R) Homo sapiens This gene encodes a G protein-coupled receptor that is closely related to the glucagon receptor and binds to glucagon-like peptide-2 (GLP2). Signalling through GLP2 stimulates intestinal growth and increases villus height in the small intestine, concomitant with increased crypt cell proliferation and decreased enterocyte apoptosis. [provided by RefSeq, Dec 2014],