

G137B Polyclonal Antibody

Cat No: HR1AP12601

For research use only

Overview

Product Name	G137B Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	GPR137B TM7SF1
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Belongs to the GPR137 family.,tissue specificity:Expressed in kidney, heart, brain and placenta.,