

## GP182 Polyclonal Antibody

Cat No: HR1AP12600

For research use only

### Overview

Product Name	GP182 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	44kDa
GeneID?Human?	GPR182 ADMR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G protein-coupled receptor 182(GPR182) Homo sapiens Adrenomedullin is a potent vasodilator peptide that exerts major effects on cardiovascular function. This gene encodes a seven-transmembrane protein that belongs to the family 1 of G-protein coupled receptors. Studies of the rat counterpart suggest that the encoded protein may function as a receptor for adrenomedullin. [provided by RefSeq, Jul 2008],