



GPR78 Polyclonal Antibody

Cat No: HR1AP12591

For research use only

Overview

Product Name	GPR78 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	39kDa
GeneID?Human?	GPR78 UNQ5925/PRO19818
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>G protein-coupled receptor 78(GPR78) Homo sapiens The protein encoded by this gene belongs to the G protein-coupled receptor family, which contain 7 transmembrane domains and transduce extracellular signals through heterotrimeric G proteins. This is an orphan receptor, which displays significant level of constitutive activity. Association analysis shows preliminary evidence for the involvement of this gene in susceptibility to bipolar affective disorder and schizophrenia. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Nov 2011].</p>