

NPBW2 Polyclonal Antibody

Cat No: HR1AP12592

For research use only

Overview

Product Name	NPBW2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	NPBWR2 GPR8
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>neuropeptides B/W receptor 2(NPBWR2) Homo sapiens The protein encoded by this gene is an integral membrane protein and G protein-coupled receptor. The encoded protein is similar in sequence to another G protein-coupled receptor (GPR7), and it is structurally similar to opioid and somatostatin receptors. This protein binds neuropeptides B and W. This gene is intronless and is expressed primarily in the frontal cortex of the brain. [provided by RefSeq, Jul 2008],</p>