

GPR61 Polyclonal Antibody

Cat No: HR1AP12588

For research use only

Overview

Product Name	GPR61 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	49kDa
GeneID?Human?	GPR61 BALGR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G protein-coupled receptor 61(GPR61) Homo sapiens This gene belongs to the G-protein coupled receptor 1 family. G protein-coupled receptors contain 7 transmembrane domains and transduce extracellular signals through heterotrimeric G proteins. The protein encoded by this gene is most closely related to biogenic amine receptors. [provided by RefSeq, Jul 2008],