

GPR26 Polyclonal Antibody

Cat No: HR1AP12580

For research use only

Overview

Product Name	GPR26 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	37kDa
GeneID?Human?	GPR26
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G protein-coupled receptor 26(GPR26) Homo sapiens This gene encodes a G protein-couple receptor protein. G-protein-coupled receptors are a large family of membrane proteins that are involved in cellular responses to environmental stimuli, neurotransmitters, and hormones. The encoded protein may play a role in neurodegenerative diseases. Epigenetic silencing of this gene has been observed in gliomas. [provided by RefSeq, Sep 2016],