

GP161 Polyclonal Antibody

Cat No: HR1AP12578

For research use only

Overview

Product Name	GP161 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	58kDa
GeneID?Human?	GPR161
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	G protein-coupled receptor 161(GPR161) Homo sapiens Upon ligand binding, G protein-coupled receptors, such as GPR161, activate cytoplasmic G proteins (see GNAS, MIM 139320), allowing the receptors to transduce extracellular signals across the plasma membrane into the cell. Phosphorylation of the receptor attenuates signaling (Matteson et al., 2008 [PubMed 18250320]).[supplied by OMIM, Aug 2008],