

CALRL Polyclonal Antibody

Cat No: HR1AP12567

For research use only

Overview

Product Name	CALRL Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	50kDa
GeneID?Human?	CALCRL CGRPR
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Receptor for calcitonin-gene-related peptide (CGRP) together with RAMP1 and receptor for adrenomedullin together with RAMP2 or RAMP3 (By similarity). The activity of this receptor is mediated by G proteins which activate adenylyl cyclase.,similarity:Belongs to the G-protein coupled receptor 2 family.,subunit:Heterodimer of CALCRL and RAMP1, RAMP2 or RAMP3.,tissue specificity:Predominantly expressed in the lung and heart.,