



BKRB2 Polyclonal Antibody

Cat No: HR1AP12565

For research use only

Overview

Product Name	BKRB2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	43kDa
GeneID?Human?	BDKRB2 BKR2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>bradykinin receptor B2(BDKRB2) Homo sapiens This gene encodes a receptor for bradykinin. The 9 aa bradykinin peptide elicits many responses including vasodilation, edema, smooth muscle spasm and pain fiber stimulation. This receptor associates with G proteins that stimulate a phosphatidylinositol-calcium second messenger system. Alternate start codons result in two isoforms of the protein. [provided by RefSeq, Jul 2008],</p>