

CS012 Polyclonal Antibody

Cat No: HR1AP12559

For research use only

Overview

Product Name	CS012 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	C19orf12
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>chromosome 19 open reading frame 12(C19orf12) Homo sapiens This gene encodes a small transmembrane protein. Mutations in this gene are a cause of neurodegeneration with brain iron accumulation-4 (NBIA4), but the specific function of the encoded protein is unknown. Alternatively spliced transcript variants encoding multiple isoforms have been observed for this gene. [provided by RefSeq, Dec 2011],</p>