

## ZCCHV Polyclonal Antibody

Cat No: HR1AP12558

For research use only

### Overview

Product Name	ZCCHV Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	99kDa
GeneID?Human?	ZC3HAV1 ZC3HDC2 PRO1677
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	zinc finger CCCH-type containing, antiviral 1(ZC3HAV1) Homo sapiens This gene encodes a CCCH-type zinc finger protein that is thought to prevent infection by retroviruses. Studies of the rat homolog indicate that the protein may primarily function to inhibit viral gene expression and induce an innate immunity to viral infection. Alternative splicing occurs at this locus and two variants, each encoding distinct isoforms, are described. [provided by RefSeq, Jul 2008],