

NIPA1 Polyclonal Antibody

Cat No: HR1AP12553

For research use only

Overview

Product Name	NIPA1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	36kDa
GeneID?Human?	NIPA1 SPG6
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>non imprinted in Prader-Willi/Angelman syndrome 1(NIPA1) Homo sapiens This gene encodes a magnesium transporter that associates with early endosomes and the cell surface in a variety of neuronal and epithelial cells. This protein may play a role in nervous system development and maintenance. Multiple transcript variants encoding different isoforms have been found for this gene. Mutations in this gene have been associated with autosomal dominant spastic paraplegia 6. [provided by RefSeq, Nov 2008],</p>