

REEP1 Polyclonal Antibody

Cat No: HR1AP12551

For research use only

Overview

Product Name	REEP1 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	22kDa
GeneID?Human?	REEP1 C2orf23
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	receptor accessory protein 1(REEP1) Homo sapiens This gene encodes a mitochondrial protein that functions to enhance the cell surface expression of odorant receptors. Mutations in this gene cause spastic paraplegia autosomal dominant type 31, a neurodegenerative disorder. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Sep 2009],