

RN185 Polyclonal Antibody

Cat No: HR1AP12546

For research use only

Overview

Product Name	RN185 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse,Rat
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	21kDa
GeneID?Human?	RNF185
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	similarity:Contains 1 RING-type zinc finger.,