

## ZFY26 Polyclonal Antibody

Cat No: HR1AP12543

For research use only

### Overview

Product Name	ZFY26 Polyclonal Antibody
Source	Rabbit
Applications	IHC-p
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	279kDa
GeneID?Human?	ZFYVE26 KIAA0321
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	zinc finger FYVE-type containing 26(ZFYVE26) Homo sapiens This gene encodes a protein which contains a FYVE zinc finger binding domain. The presence of this domain is thought to target these proteins to membrane lipids through interaction with phospholipids in the membrane. Mutations in this gene are associated with autosomal recessive spastic paraplegia-15. [provided by RefSeq, Oct 2008],