

MAVS Polyclonal Antibody

Cat No: HR1AP12540

For research use only

Overview

Product Name	MAVS Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	59kDa
GeneID?Human?	MAVS IPS1 KIAA1271 VISA
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	mitochondrial antiviral signaling protein(MAVS) Homo sapiens This gene encodes an intermediary protein necessary in the virus-triggered beta interferon signaling pathways. It is required for activation of transcription factors which regulate expression of beta interferon and contributes to antiviral immunity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Sep 2011],