

CISD2 Polyclonal Antibody

Cat No: HR1AP12539

For research use only

Overview

Product Name	CISD2 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	14kDa
GeneID?Human?	CISD2 CDGSH2 ERIS ZCD2
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	CDGSH iron sulfur domain 2(CISD2) Homo sapiens The protein encoded by this gene is a zinc finger protein that localizes to the endoplasmic reticulum. The encoded protein binds an iron/sulfur cluster and may be involved in calcium homeostasis. Defects in this gene are a cause of Wolfram syndrome 2. [provided by RefSeq, Mar 2011],