

# FGGY Polyclonal Antibody

Cat No: HR1AP12537

For research use only

## Overview

Product Name	FGGY Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	60kDa
GeneID?Human?	FGGY
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>FGGY carbohydrate kinase domain containing(FGGY) Homo sapiens This gene encodes a protein that phosphorylates carbohydrates such as ribulose, ribitol, and L-arabinitol. Genome-wide association studies in some populations have found an association between polymorphisms in this gene and sporadic amyotrophic lateral sclerosis, but studies of other populations have not been able to replicate this association. Alternative splicing results in multiple transcript variants. [provided by RefSeq, May 2013],</p>