

TMM74 Polyclonal Antibody

Cat No: HR1AP12535

For research use only

Overview

Product Name	TMM74 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	33kDa
GeneID?Human?	TMEM74
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	TMEM74 (Transmembrane Protein 74) is a Protein Coding gene. An important paralog of this gene is TMEM74B. Plays an essential role in autophagy. TMEM74-induced autophagy may involve PI3K signal transduction.