

# ANO10 Polyclonal Antibody

Cat No: HR1AP12536

For research use only

## Overview

Product Name	ANO10 Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	72kDa
GeneID?Human?	ANO10 TMEM16K
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	<p>anoctamin 10(ANO10) Homo sapiens The transmembrane protein encoded by this gene is a member of a family of calcium-activated chloride channels. Defects in this gene may be a cause of autosomal recessive spinocerebellar ataxia-10. Several transcript variants encoding different isoforms have been found for this gene.[provided by RefSeq, Mar 2011],</p>