

# F176A Polyclonal Antibody

Cat No: HR1AP12534

For research use only

## Overview

Product Name	F176A Polyclonal Antibody
Source	Rabbit
Applications	WB,ELISA
Species Reactivity	Human,Mouse
Recommended Dilutions	
Immunogen	
Species	Rabbit
Storage	-20°C/1 year
Isotype	
Clonality	
Concentration	1 mg/ml
Observed band	16kDa
GeneID?Human?	FAM176A TMEM166 SP24
Human Swiss-Prot No.	
Cellular localization	
Alternative Names	
Background	function:Acts as a regulator of programmed cell death, mediating both autophagy and apoptosis.,similarity:Belongs to the FAM176 family.,tissue specificity:Expressed in lung, kidney, liver, pancreas, placenta, but not in heart and skeletal muscle.,